

Why Your IT Team Will Love Trio

Trio is designed so that when nursing is set up for success, IT's job gets easier – not harder. Here's what that looks like for your team:

- 1 Fewer Tickets. Fewer Calls.**


When nurses log in with a PIN and the cart automatically restores their height, lighting, and access settings, the most common support requests drop. Every nurse gets a cart configured for them, on any cart in the fleet, **without IT involvement.**
- 2 Active Directory and Imprivata Integration**

With N-Sight Pro, Trio supports badge-tap login through Active Directory and Imprivata. **One tap: cart adjusts, storage unlocks, EMR session opens.** IT manages users, PINs, and access centrally – no cart-by-cart configuration.
- 3 Fleet Visibility You Don't Have to Walk the Floor For**

N-Sight Fleet Management comes standard with full-featured Trio. Monitor battery health, cart utilization, and cart location across your entire fleet in real time – no client software to install. Identify at-risk carts before they become nursing complaints.
- 4 Power That Doesn't Become Your Problem**

GoLiFe's three-tier battery system (20Ah, 40Ah, 51.2Ah) lets you match power capacity to each unit's technology draw. Swappable external batteries keep carts in circulation without downtime, and N-Sight monitors battery health fleet-wide.
- 5 One Platform. Every Department.**

Trio is a single cart architecture that IT configures for **documentation, telehealth, registration, and medication workflows – rather than managing different cart models across departments.**



See Trio in action.

Schedule a virtual demo, request a rep visit, or put a demo unit on your floor.

 www.CapsaHealthcare.com/experience-trio