

Why Your Nurses Will Love Trio

Trio is built around one idea: the cart should disappear into the workflow. Here's what your nursing team will experience when it does:

1 The Cart Remembers Them

Every nurse enters their PIN (or taps their badge with N-Sight Pro) and the cart adjusts to their saved work surface height and lighting preferences – automatically, on any cart in the fleet. No manual adjustments. No inheriting the last user's setup. The cart is theirs the moment they log in.

2 It Fits Their Bodies, Not the Other Way Around


Work surface height ranges from 33.5" to 49.5" accommodates your entire team. Memory presets mean the height is right every time without stopping to adjust. N-Stride steer assist reduces push and turning force so the cart moves with them through doorways, hallways, and tight patient rooms – not against them. Twelve-hour shifts are hard enough without the cart adding to the physical load.

3 Power That Lasts the Whole Shift

GoLiFe's scalable power system means the cart finishes the shift with them. Three battery tiers (up to 27.6 hours combined runtime at the highest capacity) cover eight-, ten-, and twelve-hour shifts without a charging detour. They don't think about batteries. They don't hunt for outlets. **They work.**

4 Storage That Matches Their Workflows

Trio's configurable storage – electronic locking and non-locking bins in multiple sizes – means the cart carries what they need for their specific unit and workflow. Fewer trips away from the patient to retrieve supplies. Secure medication access when they need it. The cart is stocked for the way they actually work.



See Trio in action.

Schedule a virtual demo, request a rep visit, or put a demo unit on your floor.

 www.CapsaHealthcare.com/experience-trio